

**August 2009**

## The SFI-Bar<sup>®</sup>

**The new ingenious stress-free bar system on implants**

The SFI-Bar System is only for Sale by Southern Implants in South Africa



- \* Tension free, excellent and stable fit.
- \* Ingenious, thanks to telescope-like connections and individual shortening.
- \* Possible to fit the SFI-Bar<sup>®</sup> 2-implant directly in the mouth.



The SFI-Bar<sup>®</sup> is the innovative bar solution for removable dentures on implants in both upper and lower jaws. The bar ensures stress-free hold of the prosthesis on the implants, thus increasing the patient's comfort. All components are prefabricated. Therefore problems associated with soldering and welding conventional bars are eliminated completely.

Moreover due to its simplicity the SFI-Bar<sup>®</sup> can be easily and efficiently adapted and placed aspiration-safe to the individual patient's mouth situation chair side. Its functional design ensures great flexibility and is indicated in most cases. With the appropriate abutments, the SFI-Bar<sup>®</sup> can be used with almost any implant system.

The SFI-Bar<sup>®</sup>2-implant and 4-implant are the standard versions. These can be upgraded with the SFI-Bar<sup>®</sup> Add-on Kit (0500 0668) to solutions for 3, 5 and 6 implants.



**International College  
of Prosthodontists**

**13th Biennial Meeting**

**PLEASE COME AND JOIN US AT THE SOUTHERN  
BOOTH AT THE ICP'S 13th BIENNIAL MEETING IN  
CAPE TOWN; SA FROM 10-13 SEPTEMBER 2009**

The ICP was founded to meet the global needs for prosthodontists and their patients, the College is an organ for international information transfer through its biennial meetings, journal, and personal communication. The international camaraderie developed through College meetings allows the specialist to share problems and solutions with colleagues who have similar concerns and goals. We are working to make the services of prosthodontists better understood by our patients and to advance the specialty globally.

For more info please visit:  
[www.southernimplants.com](http://www.southernimplants.com)  
[www.icp-org.com](http://www.icp-org.com)

CAT-0057-00 (C400)

**The Brightest Ideas in Dental Implants**

**South Africa**  
Southern Implants (Pty) Ltd.  
info@southernimplants.com

**Americas / Asia**  
Southern Implants Inc.  
info@southernimplants.us

**United Kingdom**  
Southern Implants UK  
hilary@southern-ventura.co.uk

**Australia**  
Henry Schein | Halas Dental  
info@henryschein.com.au

**New Zealand**  
Southern Implants Ltd NZ  
dkshep@xtra.co.nz

**Germany**  
Southern Implants  
info@southernimplants.de

**Greece**  
Southern Implants  
info@southernimplants.gr

**Benelux**  
ProScan bvba  
info@proscan.be

**Nordic Countries**  
Protera AB  
villy@protera.se

**Spain / Portugal**  
Contactodent  
southernimplants@sapo.pt

**Turkey**  
Ekodent  
selvasoyturk@ekodent.com

**Namibia**  
Skydancer  
mfdental@iway.na